

INSTRUCTION OF COTTON SWAB


Before use the COTTON SWAB, please read the below instruction carefully.

1. Production type/size

G1022 COTTON SWABS WITH PLAIN TUBE

2. Structure

As the below picture shows of the cotton swab structure



3. Products appearance

The tube should feel smoothly, intact, don't allow to have defects, such as sand holes, rub mark, air bubble, dirty, impurity ,etc. The tube mouth side should feel smoothly, do not have burr or flashing. The tube bottom should feel smoothly, high thickness uniformity. Transparent tube should colorless and transparent, mesopic vision can see the sample clearly and ametabolic. Swab cotton ball should be clear, no smudginess, yellowing or blackspot. Swab and tube stopper connection should match well. No gravity loose or fall off. Cotton ball of swab should wind stick tightly, No loose or fall off. Tube with swab should be sterilized, the item is EO sterile, The remaining of EO shouldn't more than 10mg/k.

4.Usage :

Gently unscrew the swab plug and remove the swab. After collecting human secretions with cotton head, the swab is inserted into the test tube and the swab plug is pressed well for laboratory use;

5. Application of scope

For clinical collection of human secretions

6. Storage transport

1、 Do not impact and extrude during transportation to avoid test tube rupture and swab rupture.

2、 The swab should be stored in a non-corrosive gas and well ventilated environment with a relative humidity of not more than 80%.

7. Expiry

A packed swab shall have a valid life of 2 years, subject to the rules of storage, transportation, storage and use.

8. Notes

- 1、 The test tube and the swab with the seal broken cannot be used.
- 2、 A swab with moldy and discolored cotton head can't be used.
- 3、 The swab may be broken during use. When inserting into human body, the swab should not be too hard to avoid breaking the skin or other harm to patients.
- 4、 Cotton swab cannot be used when it has cracks
- 5、 Products beyond the expiry date cannot be used.

9. Warning

- 1、 Sterilized products, damaged packaging should not be used.
- 2、 This product is disposable medical supplies. Once used, it must be recycled and destroyed according to regulations.

10、 Symbol explanation

Symbol	Explain	Symbol	Explain
	Can't be used second time		Can't be rained
	Easily broken, handle carefully		UP

Contact:

Zhejiang Gongdong Medical Technology Co., Ltd

ADD: No.10 Beiyuan Ave., HuangYan, Taizhou,Zhejiang, CHINA,318020

TEL: 0086-576-84082907 E-mail: stephen@chinagongdong.com

FAX: 0086-576-84050789 URL: <http://www.chinagongdong.com>